

Anti-NGFR [R05-7D3] Monoclonal Antibody

AGMB04070

Description

This Anti-NGFR [R05-7D3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04070
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	NGFR, TNFRSF16, Tumor necrosis factor receptor superfamily member 16, Gp80-LNGFR, Low affinity neurotrophin receptor p75NTR, Low-affinity nerve growth factor receptor, NGF receptor, p75 ICD, CD antigen CD271
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R05-7D3
Clonality:	Monoclonal Antibody
SwissProt ID:	P08138
Immunogen:	A synthesized peptide derived from human p75 NGF Receptor
Gene ID:	4804

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: NGFR

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 45 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IHC-P	1:100-1:200
	ICC/IF	1:20-1:50
	FC	1:50-1:100
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.