

Anti-S6K1 [R07-7V3] Monoclonal Antibody

AGMB04077

Description

This Anti-S6K1 [R07-7V3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04077
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	RPS6KB1, STK14A, Ribosomal protein S6 kinase beta-1, S6K-beta-1, S6K1, 70 kDa ribosomal protein S6 kinase 1, P70S6K1, p70-S6K 1, Ribosomal protein S6 kinase I, Serine/threonine-protein kinase 14A, p70 ribosomal S6 kinase alpha, p70 S6 kinas, KS6B1_HUMAN, p70 alpha,
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R07-7V3
Clonality:	Monoclonal Antibody
SwissProt ID:	P23443
Immunogen:	A synthesized peptide derived from human p70 S6 Kinase

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	6198
Gene Name:	RPS6KB1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 59 kDa, Observed MW: 70 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:1000-1:2000
IHC	1:100-1:200
IF	1:50-1:200
FC	1:50-1:100
IP	1:50

Protein Quantification (Optional):

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.