

Anti-AR [R06-6I9] Monoclonal Antibody

AGMB04166

Description

This Anti-AR [R06-6I9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04166
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Androgen Receptor,KD, AIS, AR8, TFM, DHTR, SBMA, HYSP1, NR3C4, SMAX1, HUMARA, Androgen Receptor
Applications:	WB IHC-P ICC/IF IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R06-6I9
Clonality:	Monoclonal Antibody
SwissProt ID:	P10275
Immunogen:	A synthesized peptide derived from human Androgen Receptor
Gene ID:	367

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: AR
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 99 kDa, Observed MW: 110 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IF	1:50-1:200
	IP	1:30

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.