

Anti-PRMT6 (2C3) [2C3-F1-G1] Monoclonal Antibody

AGMB04196

Description

This Anti-PRMT6 (2C3) [2C3-F1-G1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04196
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	HRMT1L6, FLJ10559, FLJ51477, PRMT6, ANM6_HUMAN, Chromobox protein homolog 7, FLJ10559,
Applications:	WB IHC-F IHC-P ICC/IF
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	2C3-F1-G1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q96LA8
Immunogen:	Purified recombinant human PRMT6 protein fragments expressed in E.coli.
Gene ID:	55170

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	PRMT6
Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 42 kDa, Observed MW: 42 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
	IF	1:50-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.