

Anti-IKB beta (5H3) [5H3-4C7-1H6] Monoclonal Antibody

AGMB04225

Description

This Anti-IKB beta (5H3) [5H3-4C7-1H6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04225
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	IκB-β, NFKBIB, IKB, TRIP9, NF-kappa-B inhibitor beta, NF-kappa-BIB, I-kappa-B-beta, IκB-B, IκB-beta, IkappaBbeta, Thyroid receptor-interacting protein 9, TR-interacting protein 9, TRIP-9
Applications:	IHC-P
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	5H3-4C7-1H6
Clonality:	Monoclonal Antibody
SwissProt ID:	Q15653
Immunogen:	Purified recombinant protein expressed in E.coli.
Gene ID:	4793

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: NFKBIB

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.