

# Anti-JAK1 (8D8) [8D8-2B5-5E3] Monoclonal Antibody

## AGMB04442

### Description

---

This Anti-JAK1 (8D8) [8D8-2B5-5E3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB04442
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	JTK3, AIIDE, JAK1A, JAK1B, JAK1
<b>Applications:</b>	<b>IHC-P</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	8D8-2B5-5E3
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P23458
<b>Immunogen:</b>	Synthetic peptide conjugated to KLH.
<b>Gene ID:</b>	3716
<b>Gene Name:</b>	JAK1

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse

**Isotype:** IgG1

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:100

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.