

Anti-Oct-4 (3F5) [3F5-F9-E5] Monoclonal Antibody

AGMB04505

Description

This Anti-Oct-4 (3F5) [3F5-F9-E5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04505
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Oct-4 (3F5), Oct3, Oct4, Otf3, Otf4, NF-A3, Oct-3, Oct-4, Otf-3, Otf-4, Otf3g, Oct3/4, Oct-3/4, Otf3-rs7.
Applications:	WB
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	3F5-F9-E5
Clonality:	Monoclonal Antibody
SwissProt ID:	P20263
Immunogen:	Purified recombinant mouse POU5F1/OCT4 protein fragments expressed in E.coli.
Gene ID:	18999

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: Pou5F1

Host Species: Mouse

Isotype: IgG2b

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 38 kDa, Observed MW: 45 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

	Application	Antibody Dilution Ratio
Antibody Dilution Ratio:	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.