

Anti-Influenza A Hemagglutinin H5 (3B2) [3B2-G8-C2] Monoclonal Antibody

AGMB04662

Description

This Anti-Influenza A Hemagglutinin H5 (3B2) [3B2-G8-C2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04662
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Influenza A virus H5 HA, Hemagglutinin
Applications:	ELISA
Research Area:	Microbiology
Form:	Liquid

Antibody Data

Reactivity:	Virus
Clone:	3B2-G8-C2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q77XR3
Immunogen:	Purified avian influenza virus type A (H5N1).
Gene ID:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: H5

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	ELISA	1:10000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.