

# Anti-IGJ [7A6] Monoclonal Antibody

AGMB04690

## Description

---

This Anti-IGJ [7A6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB04690
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	IGJ, Immunoglobulin J Chain, IGJ, Joining chain of multimeric IgA and IgM
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b>
<b>Research Area:</b>	Immunology
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	7A6
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P01591
<b>Immunogen:</b>	Recombinant human J chain
<b>Gene ID:</b>	3512
<b>Gene Name:</b>	JCHAIN

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse

**Isotype:** IgG2b

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** Calculated MW: 18 kDa, Observed MW: 18 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 0.09% sodium azide as a preservative.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:100-1:200
	IHC	1:200-1:500

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.