

Anti-VASH1 [14C3] Monoclonal Antibody

AGMB04703

Description

This Anti-VASH1 [14C3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04703
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	VASH1, KIAA1036, VASH, Vasohibin-1
Applications:	WB
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	14C3
Clonality:	Monoclonal Antibody
SwissProt ID:	Q7L8A9
Immunogen:	Recombinant VASH1 protein
Gene ID:	22846
Gene Name:	VASH1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Mouse
Isotype:	IgG1
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 41 kDa, Observed MW: 41 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 0.09% sodium azide as a preservative.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.