

Anti-Glycophorin A [6C8-Z1-K9] Monoclonal Antibody

AGMB04882

Description

This Anti-Glycophorin A [6C8-Z1-K9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04882
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD235a, GYPA, GPA, Glycophorin-A, MN sialoglycoprotein, PAS-2, Sialoglycoprotein alpha,
Applications:	IHC-P FC
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	6C8-Z1-K9
Clonality:	Monoclonal Antibody
SwissProt ID:	P02724
Immunogen:	Peptide
Gene ID:	2993

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: GYPA

Host Species: Mouse

Isotype: IgG1

Purification: Protein G

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:50-1:100
	FC	1:50-1:100
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	