

Anti-TPRKB [3K1-C3-K5] Monoclonal Antibody

AGMB04938

Description

This Anti-TPRKB [3K1-C3-K5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04938
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	PRPK-binding protein, CGI-121, My019
Applications:	WB ICC/IF
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	3K1-C3-K5
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9Y3C4
Immunogen:	Peptide
Gene ID:	51002
Gene Name:	TPRKB

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Mouse
Isotype:	IgG1
Purification:	Protein A
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	Calculated MW: 20 kDa, Observed MW:20 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	ICC/IF	1:50
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	