

Anti-TRAP/CD40L [1C1-T2-I7] Monoclonal Antibody

AGMB04947

Description

This Anti-TRAP/CD40L [1C1-T2-I7] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04947
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	CD40 ligand, CD40-L, T-cell antigen Gp39, TNF-related activation protein, TRAP, CD154, CD40 ligand, membrane form, CD40 ligand, soluble form, CD40LG, CD40L, TNFSF5, TRAP
Applications:	FC
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	1C1-T2-I7
Clonality:	Monoclonal Antibody
SwissProt ID:	P29965
Immunogen:	Peptide
Gene ID:	959

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: CD40LG
Host Species: Mouse
Isotype: IgG1
Purification: Protein G
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.