

Anti-TRAF2 [4D8-G8-T3] Monoclonal Antibody

AGMB04951

Description

This Anti-TRAF2 [4D8-G8-T3] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

| | |
|-----------------------|---|
| SKU: | AGMB04951 |
| Contents: | 50µl, 100µl Bradford Reagent: 1 vial (2ml) |
| Synonyms: | TRAF2, TRAP3, TNF receptor-associated factor 2, E3 ubiquitin-protein ligase TRAF2, Tumor necrosis factor type 2 receptor-associated protein 3 |
| Applications: | WB ICC/IF |
| Research Area: | Cell Biology |
| Form: | Liquid |

Antibody Data

| | |
|----------------------|---------------------|
| Reactivity: | Human |
| Clone: | 4D8-G8-T3 |
| Clonality: | Monoclonal Antibody |
| SwissProt ID: | Q12933 |
| Immunogen: | Peptide |
| Gene ID: | 7186 |

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: TRAF2
Host Species: Mouse
Isotype: IgG2a
Purification: Protein G
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 53 kDa, Observed MW:53 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|---------------------------------|--------------------|--------------------------------|
| | WB | 1:500-1:1000 |
| | ICC/IF | 1:50 |

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.