

Anti-Acetyl-p53 (Lys382) (7E1) [7E1-10D6-9B2] Monoclonal Antibody

AGMB04990

Description

This Anti-Acetyl-p53 (Lys382) (7E1) [7E1-10D6-9B2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04990
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Cellular tumor antigen p53, Cys 51 stop, Tp53, Tumor protein p53
Applications:	IHC-P
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	7E1-10D6-9B2
Clonality:	Monoclonal Antibody
SwissProt ID:	P04637
Immunogen:	Synthetic peptide conjugated to KLH.
Gene ID:	7157
Gene Name:	TP53

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Acetylated

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.