

Anti-Acetyl-p53 (Lys382) [R07-8O2] Monoclonal Antibody

AGMB04992

Description

This Anti-Acetyl-p53 (Lys382) [R07-8O2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04992
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Antigen NY-CO-13, Cellular tumor antigen p53, Phosphoprotein p53, TP53, Tumor suppressor p53
Applications:	WB IHC-P ICC/IF
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R07-8O2
Clonality:	Monoclonal Antibody
SwissProt ID:	P04637
Immunogen:	A synthesized peptide derived from human p53 (acetyl K382)
Gene ID:	7157

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	TP53
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Acetylated
Molecular Weight:	Calculated MW: 44 kDa, Observed MW: 53 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	