

Anti-Acetyl-Histone H2A (Lys9) [R06-9O2] Monoclonal Antibody

AGMB04995

Description

This Anti-Acetyl-Histone H2A (Lys9) [R06-9O2] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB04995
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	H2AK9ac, H2A, H2A1B, H2AFM, HIST1H2A, Histone H2A.2, Histone H2A/a
Applications:	WB IHC-P ICC/IF
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-9O2
Clonality:	Monoclonal Antibody
SwissProt ID:	P04908
Immunogen:	A synthesized peptide derived from human Histone H2A (acetyl K9)
Gene ID:	3012
Gene Name:	H2AC4

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Acetylated
Molecular Weight:	Calculated MW: 14 kDa, Observed MW: 14 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IHC	1:100-1:200
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	