

Anti-Phospho-IKB alpha (Ser32) [R08-2D9] Monoclonal Antibody

AGMB05134

Description

This Anti-Phospho-IKB alpha (Ser32) [R08-2D9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05134
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	NFKBIA, IKBA, MAD3, NFKBI, NF-kappa-B inhibitor alpha, I-kappa-B-alpha, Ikb-alpha, IkappaBalpha, Major histocompatibility complex enhancer-binding protein MAD3
Applications:	
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R08-2D9
Clonality:	Monoclonal Antibody
SwissProt ID:	P25963
Immunogen:	A synthetic phosphopeptide corresponding to residues surrounding Ser32 of human IKB alpha

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	4792
Gene Name:	NFKBIA
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Phosphorylated
Molecular Weight:	Calculated MW: 36 kDa, Observed MW: 36 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IP	1:20

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.