

Anti-Phospho-ERK1/2 (Tyr222/Tyr205) (10C8) [10C8-9F6-10C8]
Monoclonal Antibody
AGMB05170

Description

This Anti-Phospho-ERK1/2 (Tyr222/Tyr205) (10C8) [10C8-9F6-10C8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05170
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	MAPK3, ERK1, ERT2, ERK-1, PRKM3, P44ERK1, P44MAPK, HS44KDAP, HUMKER1A, p44-ERK1, p44-MAPK, MAPK1, ERK, p38, p40, p41, ERK2, ERT1, ERK-2, MAPK2, PRKM1, PRKM2, P42MAPK, p41mapk, p42-MAPK.
Applications:	IHC-P
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	10C8-9F6-10C8
Clonality:	Monoclonal Antibody
SwissProt ID:	P27361/P28482
Immunogen:	Synthetic peptide conjugated to KLH.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 5595/5594

Gene Name: MAPK1/MAPK3

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Phosphorylated

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.