

Anti-Phospho-IGF1 Receptor (Tyr1166) [R06-6B1] Monoclonal Antibody AGMB05200

Description

This Anti-Phospho-IGF1 Receptor (Tyr1166) [R06-6B1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05200
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	IGF1R, Insulin-like growth factor 1 receptor, Insulin-like growth factor I receptor, IGF-I receptor, CD antigen CD221
Applications:	WB IP
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R06-6B1
Clonality:	Monoclonal Antibody
SwissProt ID:	P08069
Immunogen:	A synthetic phosphopeptide corresponding to residues surrounding Tyr1166 of human INSR/IGF-1R
Gene ID:	3480

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: IGF1R
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Phosphorylated
Molecular Weight: Calculated MW: 155 kDa, Observed MW: 95 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:5000
	IP	1:10-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.