

Anti-Phospho-Vimentin (Ser39) [R38-1F-8] Monoclonal Antibody

AGMB05308

Description

This Anti-Phospho-Vimentin (Ser39) [R38-1F-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05308
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	VIM, Vimentin, CTRCT30, Epididymis luminal protein 113, FLJ36605,
Applications:	WB
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R38-1F-8
Clonality:	Monoclonal Antibody
SwissProt ID:	P08670
Immunogen:	A synthesized peptide derived from human Phospho-Vimentin (S39)
Gene ID:	7431
Gene Name:	VIM

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Phosphorylated
Molecular Weight:	Calculated MW:54kDa, Observed MW:54kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.