

Anti-Phospho-EGFR (Tyr845) [R83-2O-6] Monoclonal Antibody

AGMB05310

Description

This Anti-Phospho-EGFR (Tyr845) [R83-2O-6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05310
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	EGFR, ERBB, ERBB1, HER1, Epidermal growth factor receptor, Proto-oncogene c-ErbB-1, Receptor tyrosine-protein kinase erbB-1
Applications:	WB ICC/IF
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R83-2O-6
Clonality:	Monoclonal Antibody
SwissProt ID:	P00533
Immunogen:	A synthesized peptide derived from human Phospho-EGFR (Y845)
Gene ID:	1956

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: EGFR
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Phosphorylated
Molecular Weight: Calculated MW: 134 kDa, Observed MW: 175 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:1000-1:2000
	IF	1:50-1:200
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	