

Anti-Phospho-AS160 (Thr642) [R83-5T-8] Monoclonal Antibody

AGMB05312

Description

This Anti-Phospho-AS160 (Thr642) [R83-5T-8] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05312
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	AS160, BUB2, CDC16, NIDDM5, Tbc1d4,
Applications:	WB
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R83-5T-8
Clonality:	Monoclonal Antibody
SwissProt ID:	O60343
Immunogen:	A synthesized peptide derived from human Phospho-AS160 (T642)
Gene ID:	9882
Gene Name:	TBC1D4

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Phosphorylated
Molecular Weight:	Calculated MW:147kDa,Observed MW:180kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 10mM PBS, pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.