

Anti-SUMO2/3 [1B9-I4-D6] Monoclonal Antibody

AGMB05353

Description

This Anti-SUMO2/3 [1B9-I4-D6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05353
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	-
Applications:	WB
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	1B9-I4-D6
Clonality:	Monoclonal Antibody
SwissProt ID:	P62805
Immunogen:	Peptide
Gene ID:	121504
Gene Name:	H4C1

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species:	Mouse
Isotype:	IgG
Purification:	Affinity Chromatography
Conjugated:	Unconjugated
Modification:	Trimethylated
Molecular Weight:	Calculated MW: 18 kDa, Observed MW: 18 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in Mouse IgG in phosphate buffered saline, glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.