

Anti-IL-1RA [R05-7M2] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05488

Description

This Anti-IL-1RA [R05-7M2] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05488
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	IL-1 Receptor Antagonist Protein, DIRA, IRAP, IL1F3, IL1RA, MVCD4, IL-1RN, IL-1ra, IL-1ra3, ICIL-1RA
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-7M2
Clonality:	Monoclonal Antibody
SwissProt ID:	P18510
Immunogen:	Recombinant protein
Gene ID:	3557

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	IL1RN
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.