

Anti-MIP-1 alpha/CCL3 [R05-5W2] Monoclonal Antibody,Capture (PBS Only)

AGMB05513

Description

This Anti-MIP-1 alpha/CCL3 [R05-5W2] Monoclonal Antibody,Capture (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05513
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	CCL3, G0S19-1, MIP1A, SCYA3, C-C motif chemokine 3, G0/G1 switch regulatory protein 19-1, Macrophage inflammatory protein 1-alpha (MIP-1-alpha), PAT 464.1, SIS-beta, Small-inducible cytokine A3, Tonsillar lymphocyte LD78 alpha protein
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R05-5W2
Clonality:	Monoclonal Antibody
SwissProt ID:	P10147
Immunogen:	Recombinant protein

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID:	6348
Gene Name:	CCL3
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purified
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.