

Anti-TREM2 [R29-9E-6] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05702

Description

This Anti-TREM2 [R29-9E-6] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05702
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	TREM2, TREM-2, Trem2a, Trem2b, Trem2c, triggering receptor expressed on myeloid cells 2
Applications:	ELISA
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Mouse
Clone:	R29-9E-6
Clonality:	Monoclonal Antibody
SwissProt ID:	Q99NH8
Immunogen:	Recombinant protein
Gene ID:	83433

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: Trem2
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: SStore at 4°C.Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.