

Anti-14-3-3 sigma/SFN [R04-4T-1] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05754

Description

This Anti-14-3-3 sigma/SFN [R04-4T-1] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05754
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	SFN, HME1, 14-3-3 protein sigma, Epithelial cell marker protein 1, Stratifin
Applications:	ELISA
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R04-4T-1
Clonality:	Monoclonal Antibody
SwissProt ID:	P31947
Immunogen:	Recombinant protein
Gene ID:	2810

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: SFN
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purification
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 0.01M PBS, pH 7.2

Antibody Dilution Ratio: Optimize antibody dilution ratio according to assay requirements and sample type.

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.