

Anti-Calcitonin [R08-6T-4] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05762

Description

This Anti-Calcitonin [R08-6T-4] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05762
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	Alpha CGRP, CALC 1, CALC A, CALC1, CALCA, Calcitonin related polypeptide alpha, CGRP 1, CGRP, CGRP I, CGRP-1, CGRPI, CGRP1, CT, Katalcalcin, KC, MGC126648, alpha type CGRP, calcitonin 1, calcitonin gene related peptide 1, calcitonin gene related peptide I, calcitonin isoform CT preproprotein.
Applications:	ELISA
Research Area:	Neuroscience
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R08-6T-4
Clonality:	Monoclonal Antibody
SwissProt ID:	P01258

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Immunogen:	Recombinant protein
Gene ID:	796
Gene Name:	CALCA
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purification
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.