

Anti-VCAM1 [R03-9L-1] Monoclonal Antibody, Unconjugated, Detector (PBS Only)

AGMB05898

Description

This Anti-VCAM1 [R03-9L-1] Monoclonal Antibody, Unconjugated, Detector (PBS Only) is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB05898
Contents:	100µg Bradford Reagent: 1 vial (2ml)
Synonyms:	V-CAM 1, VCAM-1, INCAM-100, CD106
Applications:	ELISA
Research Area:	Signal Transduction
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-9L-1
Clonality:	Monoclonal Antibody
SwissProt ID:	P19320
Immunogen:	Recombinant protein
Gene ID:	7412

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name:	VCAM1
Host Species:	Rabbit
Isotype:	IgG
Purification:	Affinity Purification
Conjugated:	Unconjugated
Modification:	Unmodified
Molecular Weight:	-

Preparation & Storage

Storage:	SSStore at 4°C.Avoid freeze/thaw cycles. Store Bradford Reagent at Room Temperature for 1 Year.
Storage Buffer:	Liquid in 0.01M PBS, pH 7.2
Antibody Dilution Ratio:	Optimize antibody dilution ratio according to assay requirements and sample type.
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.