

# Anti-beta Catenin [9A4-10E6-7G5] Monoclonal Antibody

## AGMB06057

### Description

---

This Anti-beta Catenin [9A4-10E6-7G5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB06057
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	CTNNB1, CTNNB, OK/SW-cl.35, Catenin beta-1, Beta-catenin
<b>Applications:</b>	<b>IHC-P</b> <b>ELISA</b>
<b>Research Area:</b>	Signal Transduction
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	9A4-10E6-7G5
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	P35222
<b>Immunogen:</b>	Synthesized peptide derived from human Beta-catenin
<b>Gene ID:</b>	1499
<b>Gene Name:</b>	CTNNB1

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Mouse

**Isotype:** IgG2b

**Purification:** Affinity Chromatography

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

<b>Antibody Dilution Ratio:</b>	<b>Application</b>	<b>Antibody Dilution Ratio</b>
	IHC	1:200-1:400
	ELISA	1:500-1:5000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.