

Anti-CD64 [ZN74] Monoclonal Antibody

AGMB06223

Description

This Anti-CD64 [ZN74] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06223
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FCGR1A, Fc-gamma RI, FcRI, Fc-gamma RIA, FcgammaRIa, CD64, FCG1, FCGR1, IGFR1, FCRI, CD64A, FcgammaRI, CD 64, CD64 antigen,
Applications:	IHC-P
Research Area:	Immunology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	ZN74
Clonality:	Monoclonal Antibody
SwissProt ID:	P12314
Immunogen:	Full length human recombinant protein of human FCGR1A(NP_000557) produced in HEK293T cell.
Gene ID:	2209

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: FCGR1A

Host Species: Mouse

Isotype: IgG1

Purification: Affinity Chromatography

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.