

Anti-Von Willebrand Factor [ZN279] Monoclonal Antibody

AGMB06273

Description

This Anti-Von Willebrand Factor [ZN279] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06273
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	VWF, von Willebrand factor, von Willebrand antigen II, F8VWF
Applications:	IHC-P
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	ZN279
Clonality:	Monoclonal Antibody
SwissProt ID:	P04275
Immunogen:	Human recombinant protein fragment of human VWF (NP_000543) produced in SF9 cell.
Gene ID:	7450

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: VWF
Host Species: Mouse
Isotype: IgG2b
Purification: Affinity Chromatography
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Antibody	Application	Antibody Dilution Ratio
Dilution Ratio:	IHC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.