

# Anti-Cytokeratin 6 [R05-2E-5] Monoclonal Antibody

## AGMB06371

### Description

---

This Anti-Cytokeratin 6 [R05-2E-5] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

|                       |   |
|-----------------------|---|
| <b>SKU:</b>           | AGMB06371                                     |
| <b>Contents:</b>      | 50µl, 100µl<br>Bradford Reagent: 1 vial (2ml) |
| <b>Synonyms:</b>      | K6E, KRT6E, PPKNEFD, Cytokeratin 6 (KRT6)     |
| <b>Applications:</b>  | <b>IHC-P</b>                                  |
| <b>Research Area:</b> | Signal Transduction                           |
| <b>Form:</b>          | Liquid  |

### Antibody Data

---

|                      |  |
|----------------------|--|
| <b>Reactivity:</b>   | Human  |
| <b>Clone:</b>        | R05-2E-5                                       |
| <b>Clonality:</b>    | Monoclonal Antibody                            |
| <b>SwissProt ID:</b> | P48668   |
| <b>Immunogen:</b>    | A synthesized peptide derived from human KRT6A |
| <b>Gene ID:</b>      | 286887   |
| <b>Gene Name:</b>    | KRT6C  |

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Unmodified

**Molecular Weight:** -

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

| Antibody Dilution Ratio: | Application | Antibody Dilution Ratio |
|--------------------------|-------------|-------------------------|
|                          | IHC         | 1:100-1:200             |

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.