

Anti-BAP1 [R08-5M1] Monoclonal Antibody

AGMB06482

Description

This Anti-BAP1 [R08-5M1] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06482
Contents:	50µl Bradford Reagent: 1 vial (2ml)
Synonyms:	BAP 1, Bap1, BAP1_HUMAN, UCHL2, hucep-6, HU CEP-13, BAP 1, Bap1, BAP1_HUMAN,
Applications:	IHC-P
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R08-5M1
Clonality:	Monoclonal Antibody
SwissProt ID:	Q92560
Immunogen:	A synthesized peptide derived from human BAP1.
Gene ID:	8314

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: BAP-1
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC	1:100-1:200

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.