

Anti-HNF 1 beta [8M1-D3-E9] Monoclonal Antibody

AGMB06533

Description

This Anti-HNF 1 beta [8M1-D3-E9] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06533
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	FJHN, HNF 1B, HNF 2, HNF1 beta, HNF1B, HNF2, HPC11, LFB3, MODY5, TCF 2, TCF2, vHNF1
Applications:	IHC-P
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	8M1-D3-E9
Clonality:	Monoclonal Antibody
SwissProt ID:	P35680
Immunogen:	-
Gene ID:	6928

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: HNF1B
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purification
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA. Stable for 12 months from date of receipt.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	IHC-P	1:500-1:5000

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.