

Anti-MAD2L2 [R08-7E2] Monoclonal Antibody -Knockout Validated

AGMB06541

Description

This Anti-MAD2L2 [R08-7E2] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06541
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	REV7, FANCV, MAD2B, POLZ2
Applications:	WB IHC-P ICC/IF FC IP
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat
Clone:	R08-7E2
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9UI95
Immunogen:	A synthetic peptide of human Mad2L2
Gene ID:	10459
Gene Name:	MAD2L2

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Host Species: Rabbit

Isotype: IgG

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 24 kDa, Observed MW: 24 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:20-1:50
	ICC/IF	1:50-1:100
	FC	1:50-1:100
	IP	1:20-1:50

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.