

# Anti-SIRT1 [5B2] Monoclonal Antibody -Knockout Validated

## AGMB06554

### Description

---

This Anti-SIRT1 [5B2] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

### Product Information

---

<b>SKU:</b>	AGMB06554
<b>Contents:</b>	50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	SIRT1, SIR2L1, NAD-dependent protein deacetylase sirtuin-1, hSIRT1, Regulatory protein SIR2 homolog 1, SIR2-like protein 1, hSIR2
<b>Applications:</b>	<b>WB</b> <b>IHC-P</b>
<b>Research Area:</b>	Cell Biology
<b>Form:</b>	Liquid

### Antibody Data

---

<b>Reactivity:</b>	Human, Mouse, Rat
<b>Clone:</b>	5B2
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q96EB6
<b>Immunogen:</b>	Purified recombinant human SIRT1 protein fragments expressed in E.coli.
<b>Gene ID:</b>	23411

#### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene Name:** SIRT1  
**Host Species:** Mouse  
**Isotype:** IgG1  
**Purification:** Affinity Purified  
**Conjugated:** Unconjugated  
**Modification:** Unmodified  
**Molecular Weight:** Calculated MW: 82 kDa, Observed MW: 110 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
 Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in PBS containing 50% glycerol, 1% BSA and 0.02% sodium azide, pH 7.3.

	Application	Antibody Dilution Ratio
<b>Antibody Dilution Ratio:</b>	WB	1:500-1:1000
	IHC	1:50-1:100

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.