

Anti-ACAT2/Acetyl-CoA acetyltransferase [300] Monoclonal Antibody - Knockout Validated

AGMB06700

Description

This Anti-ACAT2/Acetyl-CoA acetyltransferase [300] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06700
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	Acat2, ACTL, Acetyl-CoA transferase-like protein, Cytosolic acetoacetyl-CoA thiolase, Acat2, acetoacetyl CoA thiolase, acetoacetyl Coenzyme A thiolase,
Applications:	WB IHC-P FC
Research Area:	Cardiovascular
Form:	Liquid

Antibody Data

Reactivity:	Human, Mouse, Rat, Monkey
Clone:	300
Clonality:	Monoclonal Antibody
SwissProt ID:	Q9BWD1
Immunogen:	Full length human recombinant protein of human ACAT2 produced in 293T cell.

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 39

Gene Name: ACAT2

Host Species: Mouse

Isotype: IgG2b

Purification: Affinity Purified

Conjugated: Unconjugated

Modification: Unmodified

Molecular Weight: Calculated MW: 41 kDa, Observed MW: 41 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS containing 50% glycerol, 1% BSA and 0.02% sodium azide, pH 7.3.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC	1:50-1:100
	FC	1:50-1:100

Protein Quantification (Optional): To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.