

# Anti-Phospho-CDC6 (Ser106) [R01-0F2] Monoclonal Antibody - Knockout Validated AGMB06778

## Description

---

This Anti-Phospho-CDC6 (Ser106) [R01-0F2] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

## Product Information

---

<b>SKU:</b>	AGMB06778
<b>Contents:</b>	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
<b>Synonyms:</b>	CDC6, CDC18L, Cell division control protein 6 homolog, CDC6-related protein, Cdc18-related protein, HsCdc18, p62(cdc6), HsCDC6
<b>Applications:</b>	<b>WB</b>
<b>Research Area:</b>	Epigenetics
<b>Form:</b>	Liquid

## Antibody Data

---

<b>Reactivity:</b>	Human
<b>Clone:</b>	R01-0F2
<b>Clonality:</b>	Monoclonal Antibody
<b>SwissProt ID:</b>	Q99741
<b>Immunogen:</b>	A synthetic phosphopeptide corresponding to residues surrounding Ser106 of human Cdc6

### Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

**Gene ID:** 990

**Gene Name:** CDC6

**Host Species:** Rabbit

**Isotype:** IgG

**Purification:** Affinity Purified

**Conjugated:** Unconjugated

**Modification:** Phosphorylated

**Molecular Weight:** Calculated MW: 63 kDa, Observed MW: 63 kDa

## Preparation & Storage

---

**Storage:** Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.  
Store Bradford Reagent at Room Temperature for 1 Year.

**Storage Buffer:** Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000

**Protein Quantification (Optional):** To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.