

Anti-Bcl 10 [R03-2F4] Monoclonal Antibody -Knockout Validated

AGMB06799

Description

This Anti-Bcl 10 [R03-2F4] Monoclonal Antibody -Knockout Validated is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06799
Contents:	20µl, 50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	BCL10, CIPER, CLAP, B-cell lymphoma/leukemia 10, B-cell CLL/lymphoma 10, Bcl-10, CARD-containing molecule enhancing NF-kappa-B, CARD-like apoptotic protein, hCLAP, CED-3/ICH-1 prodomain homologous E10-like regulator, CIPER, Cellular homolog
Applications:	WB IHC-P
Research Area:	Cell Biology
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	R03-2F4
Clonality:	Monoclonal Antibody
SwissProt ID:	O95999
Immunogen:	Recombinant protein of human Bcl10

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene ID: 8915
Gene Name: BCL10
Host Species: Rabbit
Isotype: IgG
Purification: Affinity Purified
Conjugated: Unconjugated
Modification: Unmodified
Molecular Weight: Calculated MW: 26 kDa, Observed MW: 26 kDa

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
 Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Antibody Dilution Ratio:

Application	Antibody Dilution Ratio
WB	1:1000-1:5000
IHC-P	1:200-1:500

Protein Quantification (Optional):

To quantify total protein levels, use the Bradford Reagent included in this kit. Visit <https://www.assaygenie.com/bradford-protein-assay-protocol/> to view the full protocol.