

Anti-Acetyl-Histone H2B (Lys11) [4J2-N5-G6] Monoclonal Antibody

AGMB06823

Description

This Anti-Acetyl-Histone H2B (Lys11) [4J2-N5-G6] Monoclonal Antibody is supplied as a kit for advanced applications. The kit includes Bradford Reagent to quantify total protein concentration for accurate sample normalization (Optional).

Product Information

SKU:	AGMB06823
Contents:	50µl, 100µl Bradford Reagent: 1 vial (2ml)
Synonyms:	H2BK11ac, Histone H2B.1 A, Histone H2B.a (H2B/a), Histone H2B.g (H2B/g), Histone H2B.h (H2B/h), Histone H2B.k (H2B/k), Histone H2B.l (H2B/l)
Applications:	WB IHC-P CHIP
Research Area:	Epigenetics
Form:	Liquid

Antibody Data

Reactivity:	Human
Clone:	4J2-N5-G6
Clonality:	Monoclonal Antibody
SwissProt ID:	P62807
Immunogen:	Peptide
Gene ID:	-

Manufacturers Statement

This final kit system is assembled and quality-released by Assay Genie Limited.

Gene Name: H2BC4

Host Species: Mouse

Isotype: IgG

Purification: Protein G

Conjugated: Unconjugated

Modification: Acetylated

Molecular Weight: -

Preparation & Storage

Storage: Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Store Bradford Reagent at Room Temperature for 1 Year.

Storage Buffer: Liquid in PBS, Glycerol and BSA

Antibody Dilution Ratio:	Application	Antibody Dilution Ratio
	WB	1:500-1:1000
	IHC-P	1:50-1:100
	ChIP	1:20
Protein Quantification (Optional):	To quantify total protein levels, use the Bradford Reagent included in this kit. Visit https://www.assaygenie.com/bradford-protein-assay-protocol/ to view the full protocol.	