

Ranibizumab Non-Neutralizing Antibody

SKU: GXDA00037 | Target: Ranibizumab | Monoclonal | Anti-Drug Antibody

Ranibizumab Non-Neutralizing Antibody (SKU: GXDA00037) is a monoclonal anti-drug antibody (ADA) supplied by Assay Genie for research into immunogenicity against Ranibizumab. It is raised in rabbit and is intended for the detection and characterisation of anti-Ranibizumab antibody responses in ELISA. This clone is characterised as non-neutralizing and binds Ranibizumab without blocking its interaction with its target, which suits it to total-ADA and bridging assay formats.

Anti-drug antibodies are a central readout in biotherapeutic development. An immune response raised against a molecule such as Ranibizumab can alter its clearance, reduce exposure and blunt clinical response, so well-characterised ADA reagents are required as assay standards, positive controls and calibrators throughout preclinical and clinical immunogenicity testing.

Specifications

Product Name	Ranibizumab Non-Neutralizing Antibody
SKU	GXDA00037
Target Drug	Ranibizumab
Antibody Type	Monoclonal
Clone ID	SAA2855
Host Species	Rabbit
Isotype	IgG
Applications	ELISA
Available Sizes	50 µg, 100 µg, 1 mg
CAS Number	347396-82-1

Properties

Format	Liquid
Buffer	0.01M PBS pH 7.4
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Guarantee	12 months from date of dispatch

Storage & Stability

Use a manual defrost freezer and avoid repeated freeze thaw cycles. Store at 2 to 8°C for frequent use. Store at -20 to -80°C for twelve months from the date of receipt.

Notes

For research use only. Not for use in diagnostic or therapeutic procedures.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
Datasheet: <https://www.assaygenie.com/content/Antibody%20Genie/GXDA/GXDA00037.pdf> | Version 1.0, 28 July 2026