

mAb-1E7 Idiotypic Nanobody (NB61)

SKU: GXDA00306 | Target: mAb-1E7 | Monoclonal | Anti-Drug Antibody

mAb-1E7 Idiotypic Nanobody (NB61) (SKU: GXDA00306) is a monoclonal anti-drug antibody (ADA) supplied by Assay Genie for research into immunogenicity against mAb-1E7. It is raised in alpaca and is intended for the detection and characterisation of anti-mAb-1E7 antibody responses in ELISA.

Anti-drug antibodies are a central readout in biotherapeutic development. An immune response raised against a molecule such as mAb-1E7 can alter its clearance, reduce exposure and blunt clinical response, so well-characterised ADA reagents are required as assay standards, positive controls and calibrators throughout preclinical and clinical immunogenicity testing.

Specifications

Product Name	mAb-1E7 Idiotypic Nanobody (NB61)
SKU	GXDA00306
Target Drug	mAb-1E7
Antibody Type	Monoclonal
Clone ID	NB61
Host Species	Alpaca
Isotype	VHH-8His-Cys-tag
Applications	ELISA
Available Sizes	50 µg, 100 µg

Properties

Format	Liquid
Buffer	0.01M PBS, pH 7.4.
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Guarantee	12 months from date of dispatch

Storage & Stability

Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.

Notes

For research use only. Not for use in diagnostic or therapeutic procedures.

Assay Genie Ltd, 25 Windsor Place, Dublin 2, D02 VY42, Ireland | +353 1 563 9720 | www.assaygenie.com | techsupport@assaygenie.com
Datasheet: <https://www.assaygenie.com/content/Antibody%20Genie/GXDA/GXDA00306.pdf> | Version 1.0, 28 July 2026