

Biotinylated Anti-bovine CCL2 Polyclonal Antibody

KPAB0012

Description

This Biotinylated Anti-bovine CCL2 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU:	KPAB0012
Concentration:	1.0 mg/mL
Conjugated:	Biotin
Applications:	ELISA, ELISpot

Antibody Data

Host:	Rabbit
Reactivity:	Bovine
Immunogen:	Recombinant Bovine CCL2
Molecular Weight:	8.8 kDa (calculated) for bovine CCL2
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Dolphin CCL2: Strong; Ovine/Caprine CCL2: Strong
Homology:	Bison bison bison (American buffalo), Bos indicus (zebu), Bos mutus (wild yak), Bos taurus (cattle)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.