

## Anti-swine IFN beta Polyclonal Antibody

KPAB0021

### Description

---

This Anti-swine IFN beta Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0021
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Swine
<b>Immunogen:</b>	Recombinant Porcine IFN beta
<b>Molecular Weight:</b>	19.5 kDa (calculated) for swine IFN- $\beta$
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Bovine IFN-B: Weak; Canine IFN-B: Weak; Equine IFN-B: Weak; Mouse IFN-B: Weak
<b>Homology:</b>	Phacochoerus africanus (Common warthog), Potamochoerus larvatus (bushpig), Potamochoerus porcus (red river hog), Sus scrofa (pig)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.