

Biotinylated Anti-bovine IL-13 Polyclonal Antibody

KPAB0023

Description

This Biotinylated Anti-bovine IL-13 Polyclonal Antibody has been validated for ELISA, ELISpot applications.

Product Information

SKU: KPAB0023
Conjugated: Biotin
Applications: ELISA, ELISpot

Antibody Data

Host: Rabbit
Reactivity: Bovine
Immunogen: Recombinant Bovine IL-13
Purification: Antigen-affinity Purification
Homology: Bos indicus (zebu), Bos taurus (cattle), Bubalus bubalis (water buffalo)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.