

Anti-rabbit IL-8 Polyclonal Antibody

KPAB0043

Description

This Anti-rabbit IL-8 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0043
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Rabbit
Immunogen:	Recombinant Rabbit IL-8
Molecular Weight:	9.1 kDa (calculated) for rabbit IL-8
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN ₃
Cross-Reactivity:	Dolphin IL-8: Weak; Guinea Pig IL-8: Strong; Human IL-8: Strong; Swine IL-8: Weak
Homology:	Oryctolagus cuniculus (rabbit)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.