

## Anti-swine CCL3L1 Polyclonal Antibody

**KPAB0083**

### Description

---

This Anti-swine CCL3L1 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

### Product Information

---

<b>SKU:</b>	KPAB0083
<b>Concentration:</b>	1.0 mg/mL
<b>Conjugated:</b>	Unconjugated
<b>Applications:</b>	ELISA, ELISpot, FC, Neutralization, WB

### Antibody Data

---

<b>Host:</b>	Rabbit
<b>Reactivity:</b>	Swine
<b>Immunogen:</b>	Recombinant Porcine CCL3L1
<b>Molecular Weight:</b>	7.8 kDa (calculated) for swine CCL3L1
<b>Purification:</b>	Antigen-affinity Purification
<b>Formulation:</b>	1X PBS buffer containing 0.09% NaN <sub>3</sub>
<b>Cross-Reactivity:</b>	Equine CCL3: Moderate
<b>Homology:</b>	Sus scrofa (pig)

## Preparation & Storage

**Storage:** Stable for up to twelve months from date of receipt at 2-8°C.

**Shipping:** Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.