

Anti-bovine IL-15 Polyclonal Antibody

KPAB0118

Description

This Anti-bovine IL-15 Polyclonal Antibody has been validated for ELISA, ELISpot, FC, Neutralization, WB applications.

Product Information

SKU:	KPAB0118
Concentration:	1.0 mg/mL
Conjugated:	Unconjugated
Applications:	ELISA, ELISpot, FC, Neutralization, WB

Antibody Data

Host:	Goat
Reactivity:	Bovine
Immunogen:	Recombinant Bovine IL-15
Molecular Weight:	13.2 kDa (calculated) for bovine IL-15
Purification:	Antigen-affinity Purification
Formulation:	1X PBS buffer containing 0.09% NaN3
Cross-Reactivity:	Equine IL-15: Weak; Guinea Pig IL-15: Weak; Swine IL-15: Weak
Homology:	Bos indicus (zebu), Bos mutus (wild yak), Bos taurus (cattle), Bubalus bubalis (water buffalo), Bubalus kerabau (carabao), Moschus berezovskii (Chinese forest musk deer)

Preparation & Storage

Storage: Stable for up to twelve months from date of receipt at 2-8°C.

Shipping: Shipped under temperature-controlled conditions using insulated packaging to maintain product stability.